



The Instrument Cleaner You Need for One Less Worry About Cross-Contamination

Proper cleaning of medical devices is the critical step in decontamination, prior to disinfection or sterilization. KIMBERLY-CLARK* Multi-Enzyme Detergent rapidly removes contaminants from reusable medical devices used throughout your facility.

The synergistic formulation of enzymes plus detergent works quickly, over a range of water temperatures. It penetrates and breaks down the most common contaminants found on your devices, including protein, fat and carbohydrate-based soils. KIMBERLY-CLARK Multi-Enzyme Detergent helps assure your devices are clean before being disinfected or sterilized, which can result in better reprocessing outcomes and reduced risk of cross-contamination from instruments.

KIMBERLY-CLARK* Multi-Enzyme Detergent provides benefits that can help make your facility's cleaning process easier and more complete:

- Effective; rapidly removes all types of soils
- Compatible; can be used with a wide range of devices
- Low-Foaming; appropriate for manual and small washer use
- Easy-To-Use; pump dispenses consistent 1 ounce dose
- User-Friendly; pleasant fragrance in an ergonomically designed package
- Process-Friendly; easily integrates into your current practices
- Environment-Friendly; biodegradable and phosphate-free

KIMBERLY-CLARK* Multi-Enzyme Detergent is just one of the clinical solutions you can depend on to meet the demands of your fast-paced world.

Whether your needs involve preventing healthcare-associated infections, surgical and digestive health solutions, or pain management, with Kimberly-Clark you'll always have one less worry.

For more information about KIMBERLY-CLARK* Multi-Enzyme Detergent, ask your Kimberly-Clark Sales Representative, authorized Kimberly-Clark Distributor Representative or visit our website at www.kchealthcare.com.



Kimberly-Clark

*Trusted Clinical Solutions**

Choosing the right detergent is just the beginning. Your Kimberly-Clark Sales Representative can work with you to help assure that you are using KIMBERLY-CLARK* Multi-Enzyme Detergent to achieve optimal efficacy.



Soak Time



Water Temperature



Dilution Ratio



Commitment to Excellence

If, for any reason, our products do not meet your expectations, please let us know your comments or suggestions for improvement. Your input will result in a concerted effort on our part to meet your requirements. Our goal is to provide quality products that completely meet your needs time after time.

For more information, please call 1-800-KCHELPS in the United States, ask your sales representative, or visit our web site at www.kchealthcare.com.

* Registered Trademark or Trademark of Kimberly-Clark Worldwide, Inc. or its affiliates. © 2006 KCWW. All rights reserved.

H8554 KLM-2116

KIMBERLY-CLARK* Multi-Enzyme Detergent

Stock #	Description	Case Information
65201	1 Gallon Bottle	4 (1 gallon) bottles w/1 dispensing pump
65205	5 Gallon Container	1 (5 gallon) container

Also available: 24 oz. KIMBERLY-CLARK* Pre-Cleaning Detergent Spray

Related Products and Services Available from Kimberly-Clark

PPE PRODUCTS

- Fluid-Resistant Facial and Eye Protection
- Impervious Cover Gowns
- Nitrile Exam Gloves
- Bouffant Caps
- Shoe Covers

ENDOSCOPY PRODUCTS

- Endoscope Cleaning Brushes
- Endoscope Cleaning Sponges

EDUCATIONAL PROGRAMS

"Coming Clean: An Essential for Reusable Medical Devices"

An educational video presentation, CE-accredited for both Nurses and CS/SPD Technicians. Ask your Kimberly-Clark Sales Representative to schedule a presentation at your facility.