

Facial Protection **Child's Face Mask**



KIMBERLY-CLARK* Child's Face Mask



Facial Protection Now Available To Fit Toddlers to Tweens! Only from Kimberly-Clark.

The Centers for Disease Control and Prevention recommends that patients with undiagnosed respiratory symptoms be asked to wear a face mask in the waiting room. The Kimberly-Clark* Child's Face Mask is designed specifically to fit small faces and provide comfortable protection for ages three and up.

- Designed specifically to fit small faces.
- Designed with lightweight materials providing a cool, comfortable fit.
- Comfortable knitted earloops.
- Formable nose wire provides a custom fit.
- Natural rubber latex free.
- Features colorful, playful designs on the outer facing.
- Wall bracket available for easy access.

The Kimberly-Clark* Child's Face Mask is recommended for use with children:

- Who are immunosuppressed.
- Who have flu-like symptoms.

From protective apparel to medical devices – everything we make, everything we do, is focused on protecting healthcare professionals and patients. This means we work every day to develop and deliver quality products and to provide responsive customer service that you can depend on to meet your protection needs.



Technical Data

Particle Filtration Efficiency > 97%
Bacterial Filtration Efficiency > 96%
Differential Pressure < 1.5

KIMBERLY-CLARK* Child's Face Masks			
Stock #	Description	Masks/Box	Boxes/Case
47127	KIMBERLY-CLARK*, Child's Face Mask	75	10

Protection. For life.





Commitment to Excellence

If, for any reason, our products do not meet your expectations, please let us know your comments or suggestions for improvement. Your input will result in a concerted effort on our part to meet your requirements. Our goal is to provide quality products that completely meet your needs time after time.

For more information, please call 1-800-KCHELPS in the United States, or visit our web site at www.kchealthcare.com.

*Registered Trademark or Trademark of Kimberly-Clark Worldwide, Inc. or its affiliates.
© 2004 KCWW. All rights reserved.

KLM-5624-1 H7396

Protection. For life.



Kimberly-Clark