

KIMBERLY-CLARK* Headwear



Kimberly-Clark offers a variety of protective headwear, with styles and protection levels for use throughout the hospital – from the cafeteria to the operating room. Bouffant caps, surgical caps and surgical hoods keep workers comfortable and well-protected, and help keep you in compliance with OSHA standards.

Bouffant Caps. Lightweight and easy to wear, Kimberly-Clark Bouffant Caps are available in both a 1-layer breathable Spunbond fabric and a 3-layer SMS fabric. The caps are designed to cover the hair in no fluid contamination or in light fluid contamination environments.

Surgical Caps. Two surgical cap styles are available to protect workers who face varying degrees of risk, depending on the procedure.

- KAYCEL* fabric caps are best suited for environments in which no fluid contamination is expected.
- Protective 3-layer fabric caps are the choice when light fluid contamination from splashing or spraying is anticipated.

Surgical Hoods. Made from a protective 3-layer fabric, Kimberly-Clark Surgical Hoods come in two models – a tie-at-the-neck style for light fluid contact, and a full-coverage hood for moderate fluid contact.

From protective apparel to surgical products – everything we do, everything we make, is focused on protecting healthcare professionals and patients. This means we work every day to develop and deliver quality products and to provide responsive customer service that you can depend on to meet your protection needs.

VALUEselect® Bouffant Cap

Coverage: Hair
 Materials: 1-layer Spunbond fabric
 Exposure Anticipated*: No potential fluid contact

Stock #	Description	Size	Items/Case	Boxes/Case
69064	VALUEselect® Bouffant Cap, Blue	Regular	500	1 (Bulk)
69067	VALUEselect® Bouffant Cap, Blue	Regular	300	3
69081	VALUEselect® Bouffant Cap, White	Regular	500	1 (Bulk)
69084	VALUEselect® Bouffant Cap, White	Regular	300	3
69071	VALUEselect® Bouffant Cap, Blue	Large	500	1 (Bulk)
69074	VALUEselect® Bouffant Cap, Blue	Large	300	3
69090	VALUEselect® Bouffant Cap, White	Large	500	1 (Bulk)
69094	VALUEselect® Bouffant Cap, White	Large	300	3

Bouffant Cap

Coverage: Hair
 Materials: 1-layer Spunbond fabric
 Exposure Anticipated*: No potential fluid contact

Stock #	Description	Size	Items/Case	Boxes/Case
69054	Bouffant Cap, Blue	Regular	500	1 (Bulk)
69221	Bouffant Cap, Blue	Regular	300	3
70535	Bouffant Cap, Cushion Band, Blue	Regular	125	4
69220	Bouffant Cap, White	Regular	300	3
69222	Bouffant Cap, Green	Regular	300	3
69053	Bouffant Cap, Blue	Large	500	1 (Bulk)
69231	Bouffant Cap, Blue	Large	300	3
69230	Bouffant Cap, White	Large	300	3
69232	Bouffant Cap, Green	Large	300	3
69004	Bouffant Cap, Blue	X-Large	300	3
60212	Bouffant Cap, Cushion Band, Blue	X-Large	250	3

Protective Bouffant Cap

Coverage: Hair
 Materials: Protective 3-layer SMS fabric
 Exposure Anticipated*: Potential for light fluid contamination from splashing or spraying

Stock #	Description	Size	Items/Case	Boxes/Case
69051	Protective Bouffant Cap, Floral	Large	240	3
69055	Protective Bouffant Cap, White	X-Large	210	3
69057	Protective Bouffant Cap, Contemporary	Large	240	3
69058	Protective Bouffant Cap, Contemporary	Regular	240	3
60400	Protective Bouffant Cap, Cushion Band, Rayon, Blue	Adjustable	75	4

Non-Elastic Tie Cap

Coverage: Hair
 Materials: Protective 3-layer SMS fabric
 Exposure Anticipated*: Potential for light fluid contamination from splashing or spraying

Stock #	Description	Size	Items/Case	Boxes/Case
69060	Non-Elastic Tie Cap, Blue	Universal	300	3

Surgical Cap

Coverage: Hair
 Materials: KAYCEL® Fabric
 Exposure Anticipated*: No potential fluid contact

Stock #	Description	Size	Items/Case	Boxes/Case
69240	Surgical Cap, Blue	Universal	300	3

Protective Surgical Cap

Coverage: Hair
 Materials: Protective 3-layer SMS fabric on sides, 1-layer fabric on top
 Exposure Anticipated*: Potential for light fluid contamination from splashing and spraying

Stock #	Description	Size	Items/Case	Boxes/Case
69520	Protective Surgical Cap, Blue	Universal	300	3

Protective Surgical Hood – Tie Neck

Coverage: Hair
 Materials: Protective 3-layer SMS fabric, 1-layer fabric on top
 Exposure Anticipated*: Potential for light fluid contamination from splashing and spraying

Stock #	Description	Size	Items/Case	Boxes/Case
69110	Protective Surgical Hood, Blue	Universal	300	3

- OSHA requires the employer to evaluate the task and type of exposure anticipated and select the appropriate protective clothing that prevents contamination of an employee's skin or clothing by blood or other potentially infectious materials.



Commitment to Excellence

If, for any reason, our products do not meet your expectations, please let us know your comments or suggestions for improvement. Your input will result in a concerted effort on our part to meet your requirements. Our goal is to provide quality products that completely meet your needs time after time.

For more information, please call 1-800-KCHELPS in the United States, or visit our web site at www.kchealthcare.com.

* Registered Trademark or Trademark of Kimberly-Clark Corp. Roswell, GA 30076 or its affiliates. © 2001 KCC. All rights reserved.

Printed on CLASSIC CREST® Text, Avon Brilliant White, 100 lb. produced by Kimberly-Clark's Neenah Paper.

KL-3638-B H5121