

## KIMBERLY-CLARK\* Drapes for General Surgical Procedures



### **For the full range of general procedures, we've got you covered.**

From laparotomy and endoscopic drapes to breast, chest and thyroid drapes, Kimberly-Clark has adaptable systems that offer maximum flexibility for surgical procedures.

### **Designed for safety and convenience.**

KIMBERLY-CLARK\* drapes are developed to protect you and your patients, with:








- Maximum rating for flame resistance for use around lasers and other electrical instruments
- Fabric that's resistant to tearing, strikethrough, and abrasion, and provides secure attachment for line and tubes.
- Resistance to strikethrough from alcohol and other prep solutions
- Lower lint than other drapes, to reduce the risk of airborne bacterial transmission
- A choice of fabric reinforcements designed either to reduce instrument slippage or control fluid run-off
- Convenient tube and line management systems, including troughs and hook & loops

### **The choice is simple.**








Now it's easier than ever to find the right KIMBERLY-CLARK drape for the job every time. Just turn the page for our one-stop selection guide.















Surgical Drapes **General Procedure Drapes**

General Surgery Drapes and Packs	Sterile Code	Non-Sterile Code	Dimensions	Fenestration	Fabric Reinforcement †	Pack Configurations	Special Application / Features
<b>Laparotomy Drapes</b>							
Laparotomy Drape 	89211	79210	76" x 120"	4" x 12"	CONTROL*		Tube Holder
Laparotomy Drape 	89221	79220 79227**	104" x 76" x 120"	4" x 12"	CONTROL*	88231; 88241; 88251; 88261; 88271	Tube Holder; Armboard Covers
Laparotomy Drape 	89228	79228	100" x 72" x 124"	4" x 12"	CONTROL* Plus	88248; 88258; 88268; 88094; 88093; 88097	Tube Holder; Armboard Covers
Laparotomy Drape w/Troughs 	89234	79234	100" x 72" x 124"	4" x 12"	CONTROL* Plus		Absorbent Troughs; Tube Holder; Armboard Covers
Laparotomy Drape w/Pouches 	89235	79235	100" x 72" x 124"	4" x 12"			Two Longitudinal Fluid Collection Pouches; Instrument Pad; Tube Holder; Armboard Covers
Laparotomy Drape w/Incise 	89282	79282	100" x 72" x 120"	8" x 13" w/ Adhesive Incise	CONTROL*		Tube Holder; Armboard Cover;
<b>Pediatric Laparotomy Drapes</b>							
Pediatric Laparotomy Drape 	89241	79240	60" x 108"	3" x 8"	CONTROL*		Tube Holders

Surgical Drapes **General Procedure Drapes**

General Surgery Drapes and Packs	Sterile Code	Non-Sterile Code	Dimensions	Fenestration	Fabric Reinforcement†	Pack Configurations	Special Application / Features
<b>Pediatric Laparotomy Drapes (Cont.)</b>							
Pediatric Laparotomy Drape II 	66886	59886	74" x 124"	4" x 6"	CONTROL* Plus		Tube Holders
Pediatric Laparotomy Drape II 	66887	59887	100" x 72" x 124"	4" x 6"	CONTROL* Plus		Armboard Covers; Tube Holders
<b>Transverse Laparotomy Drapes</b>							
Transverse Laparotomy Drape 	89231		72" x 120"	14" x 4"	CONTROL*		Tube Holders
Transverse Laparotomy Drape 	89281	79291	100" x 72" x 120"	14" x 4" transverse	CONTROL*	88291	Armboard Covers; Tube Holders
<b>Laparoscopic/Cholecystectomy Drapes</b>							
Laparoscopic/Cholecystectomy Drape 	89223	79223	102" x 76" x 120"	12.5" x 13" Abdominal	CONTROL*		Line Control Flaps; Instrument Pouch; Armboard Covers; Tube Holders
Laparoscopic/Cholecystectomy Drape 	89225	79225	102" x 76" x 120"	12.5" x 13" Abdominal	CONTROL*	88593	Line Control Flaps; Armboard Covers; Tube Holders
Laparoscopic/Lateral Offset 	89229	79229	104" x 76" x 120"	12" x 13.5" Lateral Offset		88598	Hook & Loop Line Control; Instrument Pad Reinforcement; Armboard Covers; Adjustable Fenestration Lateral Offset

General Surgery Drapes and Packs	Sterile Code	Non-Sterile Code	Dimensions	Fenestration	Fabric Reinforcement†	Pack Configurations	Special Application / Features
<b>Breast/Chest Drapes</b>							
Breast Drape 	89236	79236	100" x 76" x 124"	6" x 6"	Absorbent		Armboard Covers; Hook & Loop; Line Control
Chest Drape 	89227	79339	100" x 72" x 124"	13" x 14"	CONTROL* Plus		Instrument Pad; Armboard Covers; Tube Holders
<b>Thyroid Drapes</b>							
Thyroid Drape 	89251	79251	76" x 120"	8" x 3" Transverse	CONTROL*		Tube Holders
Thyroid Drape 	89255	79255	76" x 124"	5" x 5" Diamond	CONTROL* Plus		Instrument Pad; Tube Holders
Thyroid Drape 	89258	79258	100" x 72" x 124"	5" x 5" Diamond	CONTROL* Plus		Instrument Pad/Pouch; Armboard Covers; Tube Holders
<b>Breast/Chest Drapes</b>							
Split Drape w/ out Tape 	89321	79321	76" x 120"	4" x 40" Split	CONTROL*		Tube Holders
U-Drape 	89301	79301	76" x 120"	4" x 40" Split	CONTROL*	88311; 88321	Tube Holders

General Surgery Drapes and Packs Breast/Chest Drapes (Cont.)	Sterile Code	Non-Sterile Code	Dimensions	Fenestration	Fabric Reinforcement†	Pack Configurations	Special Application / Features
Laparoscopic/Abdominal Drape 	89233	79233	100" x 76" x 124"	12" x 13"	Absorbent	88048	Clear Pockets for Instrument Tube/Line Control; Hook & Loop Line Control; Armboard Covers
Chest U-Drape 		79335	76" x 120"	12" x 30" Split	CONTROL*	88331	Tube Holders
Split Drape w/out Tape 	89079	79079	80" x 124"	4" x 44" Split	CONTROL* Plus		Tube Holders
Bar Drape 		79370	96" x 44"		CONTROL*		Tube Holders
Bar Drape w/Arm Board Covers 		79375	106" x 30" x 62"		CONTROL*		Tube Holders

† Reinforcements include: **CONTROL\* Plus**—Maximum Absorbency; **CONTROL\***—Absorbent, non-slip for instrument control; **Absorbent Fabric**—Moderate fluid absorption; **Instrument Pad**—Non-Slip

†† Available in HANDI-BIN storage system

From medical devices to surgical apparel – everything we do, everything we make is focused on protecting healthcare professionals and patients. This means we work every day to develop and deliver quality products and to provide responsive customer service that you can depend on to meet your protection needs.



### **Commitment to Excellence**

If, for any reason, our products do not meet your expectations, please let us know your comments or suggestions for improvement. Your input will result in a concerted effort on our part to meet your requirements. Our goal is to provide quality products that completely meet your needs time after time.

**For more information, please call 1-800-KC HELPS in North America, Or visit our web site at [www.khealthcare.com](http://www.khealthcare.com).**

\* Registered Trademark or Trademark of Kimberly-Clark Worldwide, Inc., Roswell, GA 30076 or its affiliates. © 2003 KCWW. All rights reserved.

Printed on CLASSIC CREST® Text, Avon Brilliant White, 100 lb. produced by Neenah Paper Division of Kimberly-Clark.

KLM-1713 H6835

Protection. For life.

