

KIMBERLY-CLARK* Drapes for Orthopedic Procedures



For the full range of orthopedic procedures, we've got you covered.

KIMBERLY-CLARK* Orthopedic drapes offer the optimum combination of features to meet the rigorous demands of orthopedic surgery, from hand procedures to hip replacements.

KIMBERLY-CLARK* OrthoARTS* drapes provide added protection.

Made with Heavy-Duty fabric, these drapes may eliminate the need for costly drape layering. They also feature CONTROL* PLUS Fabric Reinforcement for fluid management.

Designed for safety and convenience.

KIMBERLY-CLARK* drapes are developed to protect you and your patients, with:

- Maximum rating for flame resistance for use around lasers and other electrical instruments
- Fabric that's resistant to tearing, strikethrough, and abrasion; and provides secure attachment of lines and tubes
- Much lower lint than other drapes, to reduce the risk of airborne bacterial transmission
- Fluid collection pouches
- Armboard covers for additional coverage

The choice is simple.








Now it's easier than ever to find the right KIMBERLY-CLARK drape for the job every time. Just turn the page for our one-stop selection guide.










Protection. For life.







 **Kimberly-Clark**







Surgical Drapes **Orthopedic Drapes**

Orthopedic Drapes and Packs	Sterile Code	Non-Sterile Code	Dimensions	Fenestration	Fabric Reinforcement †	Pack Configurations	Special Application / Features
General Purpose Orthopedic Drapes							
ORTHOARTS* U-Drape 	89347	79347	88" x 108"	4" x 40" split	CONTROL* PLUS	88412	Heavy-Duty Fabric; Use with ORTHOARTS* Bar drapes for maximum versatility in hip, knee or shoulder procedures;
Ortho U-Drape 		79365	88" x 108"	4" x 40" split	CONTROL*	88411	Use with our Bar drapes for versatility in hip, knee or shoulder procedures
Large Ortho U-Drape 	89311	79363	88" x 125"	6" x 54" split	CONTROL*		Generous 88" x 125" for use with larger patients. Use w/ bar drape
ORTHOARTS* Bar Drape 	89348	79348	112" x 72"		CONTROL*	88412	Heavy-Duty Fabric; Use with ORTHOARTS* U-drapes for maximum versatility in hip, knee or shoulder procedures; Tube Holders
Ortho Bar Drape 	89381	79360	100" x 60"		CONTROL*	88041	Use with U-drapes; Tube Holders
Universal Extremity Drape 	89279	79279††	90" x 124"	2.5" circle	CONTROL* PLUS	88438, 88439 88424	Conformable Fenestration; Tube Holders
Upper Extremity Orthopedic Drapes							
ORTHOARTS* Body Split Drape 	89275	79275	88" x 116"	4" x 27" Split	CONTROL* PLUS	88491	Heavy Duty Fabric; use with ORTHOARTS* Shoulder Split, or one of our Bar Drapes; Tube Holders

Surgical Drapes **Orthopedic Drapes**

Orthopedic Drapes and Packs	Sterile Code	Non-Sterile Code	Dimensions	Fenestration	Fabric Reinforcement †	Pack Configurations	Special Application / Features
Upper Extremity Orthopedic Drapes (Cont.)							
ORTHOARTS* Shoulder Split Drape 	89274	79274	112" x 60"	6" x 8" dual splits	CONTROL* PLUS	88491	Heavy Duty fabric; Use with the ORTHOARTS* Body Split drape; Features dual splits for right or left shoulder procedures; Tube Holders
ORTHOARTS* Shoulder Arthroscopy Drape 	89492	79492	100" x 60"	6" x 3" x 6"		88492	Heavy Duty fabric; Attached fluid pouch; Use with ORTHOARTS* Body Split Drape
Beach Chair Shoulder Drape w/ Pouch 	89066	79049	172" x 100"		CONTROL* PLUS		Attachable fluid pouch with suction port; Complete one-drape system; Hook & loop line controls; Pouch attachment zones
Beach Chair Shoulder Drape 	89063		172" x 100"		CONTROL* PLUS	88108	Complete one-drape system; Attachable pouch available; Hook & loop line controls
Beach Chair Shoulder Drape II 	89069	79069	160" x 102"	2.5" Conformable Tear Drop	Instrument pad		One-piece drape with attached fluid-collection pouch; Hook & loop tube holders
Upper Extremity Drape 	89271	79273	76" x 120"	1.25" conformable Circular	CONTROL*	88421	Conformable fenestration for a tight seal; Tube Holders
ORTHOARTS* Hand Drape 	89272	79272	76" x 150"	1.25" conformable Circular	CONTROL* PLUS		Heavy-Duty Fabric; Include 62" x 58" gusseted hand table cover; Tube holders

Orthopedic Drapes and Packs	Sterile Code	Non-Sterile Code	Dimensions	Fenestration	Fabric Reinforcement	Pack Configurations	Special Application / Features
Upper Extremity Orthopedic Drapes (Cont.)							
Hand Drape 	77450	76379	114" x 142"	1.25" conformable circular	CONTROL* Plus		Includes 62" x 58" hand table cover; Tube holders
Hip Orthopedic Drapes							
ORTHOARTS* Hip Drape w/Armboard Covers 	89346	79346	114" x 88" X 126"	8.25" X 6.75" conformable teardrop	CONTROL* PLUS	88453	Heavy-Duty fabric; Armboard covers; Conformable fenestration; Tube holders
ORTHOARTS* Hip Drape 	89345	79345	114" x 88" X 126"	8.25" X 6.75" conformable teardrop	CONTROL* PLUS	88452	Heavy-Duty fabric; Leg pockets on both sides; Armboard covers; Tube holders
Hip Drape 	89341	79341	116" x 88" X 131"	Conformable teardrop	CONTROL*	88451 88441	Armboard covers; Conformable fenestration; Tube holders
Hip Drape w/ Clear Leg Pockets 	89349	79349	116" x 136"	Conformable teardrop	CONTROL* PLUS		Clear leg pockets; Tube holders; Armboard covers
Lower Extremity Orthopedic Drapes							
ORTHOARTS* Lower Extremity Drape 	89278	79278	114" x 88" X 126"	2.5" conformable circular	CONTROL* PLUS	88471	Conformable circular fenestration; Heavy-Duty fabric; Tube holders; Armboard cover

Orthopedic Drapes and Packs	Sterile Code	Non-Sterile Code	Dimensions	Fenestration	Fabric Reinforcement †	Pack Configurations	Special Application / Features
Lower Extremity Orthopedic Drapes (Cont.)							
Lower Extremity Drape w/ Armboard Covers 	89016	79017	108" x 80" X 124"	2.5" conformable circular	CONTROL* PLUS		Conformable circular fenestration; Tube holders; Armboard covers
Lower Extremity Drape 	89276	79276	88" x 131"	2.5" conformable circular	CONTROL*	88431, 88461, 88474	Conformable circular fenestration; Tube holders;
ORTHOARTS* Knee Arthroscopy Drape 	89267	79267††	114" x 88" x 106"	2.5" conformable circular	CONTROL* PLUS	88481	Heavy-Duty fabric; Fluid collection pouch with reflux shields to prevent backflow; Suction port; Armboard covers; Tube holders
Knee Arthroscopy Drape 	89265	79265	90" x 124"	2.5" conformable		88434, 88488	Fluid collection pouch; Tube holders; Conformable circular fenestration
Bilateral Limb Drape 	89291	79381	76" x 120"	Two 2.5" circular	CONTROL*		Two, 2.5" circular fenestrations; Tube holders
Impervious Split Drape 	89331	79330	60" x 70"	4" x 21" Split		88411; 88412; 88421; 88108; 88451; 88461; 88438; 88439; 88471; 88488; 88481; 88431	Impervious Also available in Large Impervious Split #79333 non-sterile (60" x 81") Pack 88452

† Reinforcements include: **CONTROL* Plus**—Maximum Absorbency; **CONTROL***—Absorbent, non-slip for instrument control; **Absorbent Fabric**—Moderate fluid absorption; **Instrument Pad**—Non-Slip

†† Non-Sterile in HANDI-BIN storage system

From medical devices to surgical apparel – everything we do, everything we make is focused on protecting healthcare professionals and patients. This means we work every day to develop and deliver quality products and to provide responsive customer service that you can depend on to meet your protection needs.



Commitment to Excellence

If, for any reason, our products do not meet your expectations, please let us know your comments or suggestions for improvement. Your input will result in a concerted effort on our part to meet your requirements. Our goal is to provide quality products that completely meet your needs time after time.

For more information, please call 1-800-KC HELPS in North America, Or visit our web site at www.khealthcare.com.

* Registered Trademark or Trademark of Kimberly-Clark Worldwide, Inc., Roswell, GA 30076 or its affiliates. © 2003 KCWW. All rights reserved.

Printed on CLASSIC CREST® Text, Avon Brilliant White, 100 lb. produced by Neenah Paper Division of Kimberly-Clark.

KLM-1723 H6835

Protection. For life.

