



KIMBERLY-CLARK* FLUIDSHIELD* PFR95* N95 Respirator Singles



NEW Individually Wrapped Respirators

Healthcare can be a very fast and fluid environment.

From the ER to EMS and patient care, healthcare workers need to exercise caution in the face of respiratory threats. NIOSH-approved N95 respirators are the go-to choice for airborne protection or when the actual respiratory risks are unknown.

The risk of infection with tuberculosis and other airborne diseases, such as SARS, smallpox, measles or avian influenza, is an ongoing threat to healthcare workers. The healthcare environment contains many potential respiratory hazards in the form of infectious agents from patients who release droplet nuclei into the air via coughing or other respiratory function. In addition, there is danger of infection from surgical procedures where these agents may be aerosolized.

KIMBERLY-CLARK* FLUIDSHIELD* PFR95* N95 Respirator Singles provide fluid-resistant, respiratory protection you can count on, right at your fingertips. These convenient, individually wrapped respirators are easy to carry, so you're always prepared to help yourself to respiratory protection in any setting.

- **DUCKBILL* Face Mask: Unique shape provides large breathing chamber**
- **Fold-flat packaging is easy to carry**
- **Fluid-resistant at 160 mmHg – ASTM 1862**
- **Sizes designed to provide effective and comfortable fit***
- **Ideal for use in the ER, EMS or anywhere the respiratory threat is unknown**




Trusted Clinical Solutions*

KIMBERLY-CLARK* FLUIDSHIELD* PFR95* N95 Respirator Singles

Reduce Respiratory Risk Wherever You Are:

FLUIDSHIELD* PFR95* N95 Respirator Singles are NIOSH-approved, N95 Respirators designed to provide an effective facial fit, as required for proper respiratory protection against airborne pathogens.

Code	Description	Size	Items/Box	Boxes/Case
46728	FLUIDSHIELD* PFR95* N95 Respirator Singles (IWR)	Reg.	35	6 boxes/case
46828	FLUIDSHIELD* PFR95* N95 Respirator Singles (IWR)	Sm.	35	6 boxes/case



Kimberly-Clark also has a **complete** family of nitrile exam gloves that provide protection, performance and comfort in a non-latex glove.



KC500 PURPLE NITRILE* Exam Gloves are widely used in healthcare facilities nationwide to protect clinicians from bacteria, viruses, and chemicals during procedures where the risk of fluid exposure is moderate to high. Available in 9.5" and 12" (XTRA*) lengths, these durable gloves provide maximum coverage and protection, especially when dealing with unknown risks.

Code	Description	Size	Items/Box	Boxes/Case
55080	KC500 PURPLE NITRILE* Gloves	X-Small	100 eaches/box	10 boxes/case
55081	KC500 PURPLE NITRILE* Gloves	Small	100 eaches/box	10 boxes/case
55082	KC500 PURPLE NITRILE* Gloves	Medium	100 eaches/box	10 boxes/case
55083	KC500 PURPLE NITRILE* Gloves	Large	100 eaches/box	10 boxes/case
55084	KC500 PURPLE NITRILE* Gloves	X-Large	90 eaches/box	10 boxes/case
55090	KC500 PURPLE NITRILE XTRA* Gloves	X-Small	50 eaches/box	10 boxes/case
50601	KC500 PURPLE NITRILE XTRA* Gloves	Small	50 eaches/box	10 boxes/case
50602	KC500 PURPLE NITRILE XTRA* Gloves	Medium	50 eaches/box	10 boxes/case

KC300 STERLING* NITRILE Exam Gloves provide excellent fit and feel compared to latex, excellent comfort and tactile sensitivity, and an excellent durable barrier. If you want to standardize on one exam glove that can fit virtually all your needs, facility-wide, STERLING* NITRILE is the glove to choose. Available in 9.5" or 12" (XTRA*) lengths.

Code	Description	Size	Items/Box	Boxes/Case
55085	KC300 STERLING* NITRILE Gloves	X-Small	150 eaches/box	10 boxes/case
55086	KC300 STERLING* NITRILE Gloves	Small	150 eaches/box	10 boxes/case
55087	KC300 STERLING* NITRILE Gloves	Medium	150 eaches/box	10 boxes/case
55088	KC300 STERLING* NITRILE Gloves	Large	150 eaches/box	10 boxes/case
55089	KC300 STERLING* NITRILE Gloves	X-Large	140 eaches/box	10 boxes/case
53137	KC300 STERLING* NITRILE XTRA* Gloves	X-Small	100 eaches/box	10 boxes/case
53138	KC300 STERLING* NITRILE XTRA* Gloves	Small	100 eaches/box	10 boxes/case
53139	KC300 STERLING* NITRILE XTRA* Gloves	Medium	100 eaches/box	10 boxes/case
53140	KC300 STERLING* NITRILE XTRA* Gloves	Large	100 eaches/box	10 boxes/case
53141	KC300 STERLING* NITRILE XTRA* Gloves	X-Large	100 eaches/box	10 boxes/case

The KIMBERLY-CLARK ADVANTAGE†

KNOWLEDGE NETWORK* Accredited Education
Ongoing Customer Support
Expert Sales Force
Tools & Best Practices
Clinical Research
Commitment to Excellence Guarantee

Infection prevention website:

www.HAIwatch.com



For more information, please call your sales representative, or visit our web site at www.kchealthcare.com.

*Registered Trademark or Trademark of Kimberly-Clark Worldwide, Inc. The COLOR ORANGE and COLOR PURPLE and COLOR GREY are Registered Trademarks of KCWW. ©2009 KCWW.



Trusted Clinical Solutions®

† Subject to individual Quantitative Fit Testing